


Anti-Tfe3 Polyclonal Antibody For Wb, Elisa - P19532

Item Number: CM0029843



Introduction

Rabbit polyclonal antibody against human TFE3 (Transcription Factor E3). Validated for WB and ELISA. Reacts with human, mouse, and rat TFE3. Ideal for lysosomal biogenesis and mTOR signaling research.

[Learn More](#)

Parameter	Value
Product Name	TFE3 Rabbit pAb
Remarks/Alias	TFEA; RCCP2; RCCX1; MRXSPF; bHLHe33; TFE3
Species	Human
GeneID (Human)	7030
GeneID	7030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TFE3 (NP_006512.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19532
Protein Weight	62kDa
Shipping	Ice bag