

Anti-Tfeb Polyclonal Antibody For Wb, Ip, Elisa - Q9R210

Item Number: CM0029205

Introduction

Rabbit polyclonal antibody targeting TFEB (Transcription factor EB), validated for Western blot, IP, and ELISA. Detects human and mouse TFEB, a master regulator of autophagy and lysosomal biogenesis. Ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	TFEB Rabbit pAb
Remarks/Alias	Tcfeb; bHLHe35
Species	Mouse
GeneID (Human)	1737
GeneID	21425
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 390-476 of mouse TFEB (NP_001155194.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9R210
Protein Weight	53 kDa
Shipping	Ice bag