


Anti-Tfeb Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P19484

Item Number: CM0021912



Introduction

Rabbit polyclonal antibody against human TFEB, validated for WB, IHC-P, IF/ICC, ELISA.

Recognizes human, mouse, rat TFEB. Essential tool for autophagy, lysosomal biogenesis, and transcription research.

[Learn More](#)

Parameter	Value
Product Name	TFEB Rabbit pAb
Remarks/Alias	TCFEB; BHLHE35; ALPHATFEB; TFEB
Species	Human
Gene ID (Human)	7942
Gene ID	7942
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 391-461 of human TFEB (NP_009093.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19484
Protein Weight	52kDa
Shipping	Ice bag