


Anti-Tff1/Ps2 Monoclonal Antibody For Wb - P04155

Item Number: CM0012796



Introduction

Anti-TFF1/pS2 monoclonal antibody for Western blot, immunoprecipitation, immunofluorescence, and ELISA. Recognizes human, mouse, rat trefoil factor 1. Target weight 9 kDa. Ideal for gastrointestinal and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TFF1/pS2 Rabbit PolymAb®
Remarks/Alias	pS2; BCEI; HPS2; HP1.A; pNR-2; D21S21
Species	Human
Gene ID (Human)	7031
Gene ID	7031
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04155
Protein Weight	9 kDa
Shipping	Ice bag