


Anti-Tff1/Ps2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04155

Item Number: CM0024622



Introduction

Rabbit polyclonal antibody against human TFF1/pS2 (trefoil factor 1). Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Suitable for gastrointestinal and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TFF1/pS2 Rabbit pAb
Remarks/Alias	pS2; BCEI; HPS2; HP1.A; pNR-2; D21S21; TFF1/pS2
Species	Human
Gene ID	7031
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 25-84 of human TFF1/pS2 (NP_003216.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04155
Protein Weight	9kDa
Shipping	Ice bag