

Anti-Tff2 Polyclonal Antibody For Wb, Elisa - Q03403

Item Number: CM0025750

Introduction

Rabbit polyclonal antibody specific for human Trefoil factor 2 (TFF2). Validated for WB and ELISA applications. TFF2 is a secreted protein that modulates gastric acid secretion and mucus stability, making it a key target for gastrointestinal disease research.

[Learn More](#)

Parameter	Value
Product Name	TFF2 Rabbit pAb
Remarks/Alias	SP; SML1; TFF2
Species	Human
Gene ID	7032
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-129 of human TFF2 (NP_005414.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q03403
Protein Weight	14kDa
Shipping	Ice bag