


Anti-Tff3 Monoclonal Antibody For Ihc-P, Elisa - Q07654

Item Number: CM0010176



Introduction

This rabbit monoclonal antibody targets human TFF3 (intestinal trefoil factor), validated for IHC-P and ELISA applications. Ideal for investigating mucosal repair and epithelial cell mobility in intestinal and uterine tissues.

[Learn More](#)

Specification	Details
Product Name	TFF3 Rabbit mAb
Remarks/Alias	ITF; P1B; TFI; TFF3
Species	Human
GeneID (Human)	7033
GeneID	7033
Immunogen	Synthetic peptide corresponding to amino acids 1-80 of human TFF3 (Q07654)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q07654
Protein Weight	9kDa
Shipping	Ice bag