


Anti-Tff3 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q07654

Item Number: CM0019824



Introduction

High-quality TFF3 Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence, and ELISA. Reacts with human TFF3, a key mucosal repair factor. Ideal for gastrointestinal and inflammatory bowel disease research.

[Learn More](#)

Parameter	Value
Product Name	TFF3 Rabbit pAb
Remarks/Alias	ITF; P1B; TFI; TFF3
Species	Human
Gene ID (Human)	7033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-94 of human TFF3 (NP_003217.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	Q07654
Protein Weight	9kDa
Shipping	Ice bag