

Anti-Tfpi Rabbit Monoclonal Antibody For Wb, Elisa - P10646

Item Number: CM0010722



Introduction

Rabbit monoclonal antibody targeting human TFPI (Tissue factor pathway inhibitor), validated for WB and ELISA. Essential for coagulation research and diagnostic assay development. - P10646

[Learn More](#)

Parameter	Value
Product Name	TFPI Rabbit mAb
Remarks/Alias	EPI; TFI; LACI; TFPI1; TFPI
Species	Human
GeneID (Human)	7035
GeneID	7035
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-304 of human TFPI (P10646).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10646
Protein Weight	35kDa
Shipping	Ice bag