

Anti-Tfpi2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P48307

Item Number: CM0025847



Introduction

Rabbit polyclonal antibody against human TFPI2 (PP5, TFPI-2), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Useful for studies on coagulation, matrix remodeling, and serine protease inhibition.

[Learn More](#)

Parameter	Value
Product Name	TFPI2 Rabbit pAb
Remarks/Alias	PP5, REF1, TFPI-2
Species	Human
Gene ID	7980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-235 of human TFPI2 (NP_006519.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48307
Protein Weight	27kDa
Shipping	Ice bag