


Anti-Tgf Beta 1 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01137

Item Number: CM0010908



Introduction

Rabbit monoclonal antibody against human TGF beta 1, validated for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat. Suitable for studying TGF-beta signaling, development, and immune function.

[Learn More](#)

Parameter	Value
Product Name	TGF beta 1 Rabbit mAb
Remarks/Alias	CED; LAP; DPD1; TGFB; IBDIMDE; TGFBeta; TGF-beta1; CAEND1
Species	Human
GeneID (Human)	7040
GeneID	7040
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 279-390 of human TGF beta 1 (NP_000651.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01137
Protein Weight	44 kDa
Shipping	Ice bag