


Anti-Tgf-Beta-1 Polyclonal Antibody For Wb, Elisa - P04202

Item Number: CM0023310



Introduction

Rabbit polyclonal antibody targeting mouse TGF-beta-1. Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for research in immune function, development, and cell differentiation.

[Learn More](#)

Parameter	Specification
Product Name	Tgfb1 Rabbit pAb
Remarks/Alias	Tgfb; Tgfb-1; TGFbeta1; TGF-beta1; Tgfb1
Species	Mouse
Gene ID (Mouse)	21803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-390 of mouse Tgfb1 (NP_035707.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04202
Protein Weight	44kDa
Shipping	Ice bag