

Anti-Tgf Beta 1 Polyclonal Antibody For Wb, If/Icc, Elisa - P01137

Item Number: CM0021655



Introduction

Rabbit polyclonal antibody against human TGF beta 1, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects proprotein and mature forms. Ideal for studying TGF-beta signaling, immune regulation, and EMT.

[Learn More](#)

Parameter	Value
Product Name	TGF beta 1 Rabbit pAb
Remarks/Alias	CED; LAP; DPD1; TGFB; IBDIMDE; TGFbeta; TGF-beta1; TGF beta 1
Species	Human
GeneID (Human)	7040
GeneID	7040
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 291-390 of human TGF beta 1 (NP_000651.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	P01137
Protein Weight	44kDa
Shipping	Ice bag