


Anti-Tgf Beta Receptor I (Tgfbr1) Rabbit Monoclonal Antibody For Wb And Elisa - P36897

Item Number: CM0007177



Introduction

High-quality rabbit monoclonal antibody against human TGF-beta receptor type-1 (TGFBR1/ALK5), validated for Western blot and ELISA. Ideal for studying TGF-beta signaling, cancer biology, and fibrosis research.

[Learn More](#)

Parameter	Value
Product Name	TGF beta Receptor I (TGFBR1) Rabbit mAb
Remarks/Alias	AAT5; ALK5; ESS1; LDS1; MSSE; SKR4; TBRI; ALK-5; LDS1A; LDS2A; TBR-i; TGFR-1; ACVRLK4; tbetaR-I; TGF beta Receptor I (TGFBR1)
Species	Human
GeneID (human)	7046
GeneID	7046
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human TGF beta Receptor I (TGFBR1) (NP_004603.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P36897
Protein Weight	56kDa
Shipping	Ice bag