


Anti-Tgf Beta Receptor I (Tgfbr1) Rabbit Polyclonal Antibody For Wb And Elisa - P36897

Item Number: CM0016349



Introduction

Rabbit polyclonal antibody targeting human TGFBR1, validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for TGF-beta signaling, cell proliferation, epithelial-mesenchymal transition, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TGF beta Receptor I (TGFBR1) Rabbit pAb
Remarks/Alias	AAT5; ALK5; ESS1; LDS1; MSSE; SKR4; TBRI; ALK-5; LDS1A; LDS2A; TBR-i; TGFR-1; ACVRLK4; tbetaR-I; TGF beta Receptor I (TGFBR1)
Species	Human
GeneID (human)	7046
GeneID	7046
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TGF beta Receptor I (TGFBR1) (NP_004603.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36897
Protein Weight	56 kDa
Shipping	Ice bag