


Anti-Tgf Beta Receptor Ii Monoclonal Antibody For Wb - P37173

Item Number: CM0003567



Introduction

Rabbit monoclonal antibody targeting human TGF beta Receptor II, validated for Western blot and ELISA. Ideal for TGF-beta signaling pathway research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	TGF beta Receptor II Rabbit mAb
Remarks/Alias	AAT3; FAA3; LDS2; MFS2; RIIC; LDS1B; LDS2B; TAAD2; TBRII; TBR-ii; TGFR-2; tbetaR-II; TGFbeta-RII; TGF beta Receptor II (TGFB2)
Species	Human
Gene ID	7048
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 23-130 of human TGF beta Receptor II (TGFB2) (P37173))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P37173
Protein Weight	65kDa/67kDa/9kDa
Shipping	Ice bag