


Anti-Tgf Beta Receptor Ii (Tgfbbr2) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P37173

Item Number: CM0014594



Introduction

High-quality rabbit polyclonal antibody targeting human TGF-beta receptor type-2 (TGFB2, P37173). Validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for TGF-beta signaling research in oncology and fibrosis.

[Learn More](#)

Parameter	Value
Product Name	TGF beta Receptor II Rabbit pAb
Remarks/Alias	AAT3; FAA3; LDS2; MFS2; RIIC; LDS1B; LDS2B; TAAD2; TBR11; TBR-ii; TGFR-2; tbetaR-II; TGFbeta-R11; TGF beta Receptor II (TGFB2)
Species	Human
GeneID (human)	7048
GeneID	7048
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-280 of human TGF beta Receptor II (TGFB2) (NP_003233.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37173
Protein Weight	65kDa
Shipping	Ice bag