

Anti-Tgfa Rabbit Polyclonal Antibody For Wb, Elisa - P01135

Item Number: CM0013185



Introduction

Rabbit polyclonal antibody against human TGFA, validated for WB and ELISA. Recognizes human and mouse TGFA; suitable for studies of EGFR signaling and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	TGFA Rabbit pAb
Remarks/Alias	TFGA; TGFA
Species	Human
Gene ID	7039
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-98 of human TGFA (NP_001093161.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01135
Protein Weight	17kDa
Shipping	Ice bag