

Anti-Tgfb1l1 Polyclonal Antibody For Wb, Elisa - 043294

Item Number: CM0029691

Introduction

High-quality rabbit polyclonal antibody targeting TGFB1I1 (HIC5/ARA55) for Western Blot and ELISA applications. Validated in human, mouse, and rat samples; detects a ~50kDa protein. Ideal for cell signaling and focal adhesion research.

[Learn More](#)

Parameter	Specification
Product Name	TGFB1I1 Rabbit pAb
Remarks/Alias	HIC5, ARA55, HIC-5, TSC-5, TGFB1I1
Species	Human
GeneID (Human)	7041
GeneID	7041
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-220 of human TGFB1I1 (NP_001035919.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43294
Protein Weight	50kDa
Shipping	Ice bag