

Anti-Tgfb2 Rabbit Polyclonal Antibody For Wb And Elisa - P61812

Item Number: CM0030601



Introduction

TGFB2 Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, rat. Recombinant immunogen targeting amino acids 20-160. Ideal for TGF-beta-2 signaling and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	TGFB2 Rabbit pAb
Remarks/Alias	LDS4; G-TSF; TGF-beta2; TGFB2
Species	Human
GeneID (human)	7042
GeneID	7042
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-160 of human TGFB2 (NP_001129071.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61812
Protein Weight	48kDa
Shipping	Ice bag