

Anti-Tgfb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P61812

Item Number: CM0020048



Introduction

High-quality rabbit polyclonal antibody against human TGFB2, validated for WB, IHC-P, and ELISA. Targets the TGF-beta2 proprotein, a key regulator in angiogenesis and heart development, ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	TGFB2 Rabbit pAb
Remarks/Alias	LDS4; G-TSF; TGF-beta2; TGFB2
Species	Human
Gene ID (Human)	7042
Gene ID	7042
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-160 of human TGFB2 (NP_001129071.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P61812
Protein Weight	48kDa
Shipping	Ice bag