


Anti-Tgfb3 Polyclonal Antibody For Wb, Ihc-P, Elisa - P10600

Item Number: CM0026578



Introduction

Rabbit polyclonal antibody against TGFB3, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat TGFB3. Ideal for studying TGF- β signaling, embryogenesis, and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	TGFB3 Rabbit pAb
Remarks/Alias	ARVD; LDS5; RNHF; ARVD1; TGF-beta3; TGFB3
Species	Human
Gene ID	7043
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-300 of human TGFB3 (NP_003230.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10600
Protein Weight	47kDa
Shipping	Ice bag