


Anti-Tgfb1 (Bigh3) Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15582

Item Number: CM0006405



Introduction

Rabbit monoclonal antibody targeting TGFBI (transforming growth factor-beta-induced protein ig-h3), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for cell adhesion, corneal dystrophy, and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	TGF beta induced (TGFBI) Rabbit mAb
Remarks/Alias	CSD; CDB1; CDG2; CSD1; CSD2; CSD3; EBMD; LCD1; BIGH3; CDGG1; TGF beta induced (TGFBI)
Species	Human
GeneID (human)	7045
GeneID	7045
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 550-650 of human TGFBI (NP_000349.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15582
Protein Weight	74KDa
Shipping	Ice bag