


Anti-Tgfbf Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15582

Item Number: CM0029759



Introduction

High-quality rabbit polyclonal antibody against TGFBI (BIGH3) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat TGFBI. Useful for studying cell adhesion, extracellular matrix, and corneal dystrophies.

[Learn More](#)

Parameter	Description
Product Name	TGF beta induced (TGFBI) Rabbit pAb
Remarks/Alias	CSD; CDB1; CDG2; CSD1; CSD2; CSD3; EBMD; LCD1; BIGH3; CDGG1; TGF beta induced (TGFBI)
Species	Human
Gene ID (Human)	7045
Gene ID	7045
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-650 of human TGF beta induced (TGFBI) (NP_000349.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15582
Protein Weight	75kDa
Shipping	Ice bag