


Anti-Tgfb1 (Tgf Beta Induced) Rabbit Monoclonal Antibody For Wb, Elisa - Q15582

Item Number: CM0004300



Introduction

High-quality rabbit monoclonal antibody targeting human TGFBI (TGF beta induced) for WB and ELISA. Cross-reacts with mouse and rat. Key for cell adhesion and corneal dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	TGF beta induced (TGFBI) Rabbit mAb
Remarks/Alias	CSD; CDB1; CDG2; CSD1; CSD2; CSD3; EBMD; LCD1; BIGH3; CDGG1; TGF beta induced (TGFBI)
Species	Human
Gene ID (Human)	7045
Gene ID	7045
Immunogen	A synthetic peptide corresponding to amino acids 550-650 of human TGF beta induced (TGFBI) (NP_000349.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15582
Protein Weight	75 kDa
Shipping	Ice bag