

Anti-Tgfb1 (Tgf Beta Induced) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15582

Item Number: CM0004943



Introduction

High-quality anti-TGFBI rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat TGFBI. Ideal for cell adhesion, corneal dystrophy, and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	TGF beta induced (TGFBI) Rabbit mAb
Remarks/Alias	CSD, CDB1, CDG2, CSD1, CSD2, CSD3, EBMD, LCD1, BIGH3, CDGG1, TGF beta induced (TGFBI)
Species	Human
Gene ID (Human)	7045
Gene ID	7045
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15582
Protein Weight	75 kDa
Shipping	Ice bag