

Anti-Tgif1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15583

Item Number: CM0019902



Introduction

Rabbit polyclonal antibody targeting human TGIF1 (SWISS: Q15583). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Useful for transcription regulation and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	TGIF1 Rabbit pAb
Remarks/Alias	HPE4; TGIF; TGIF1
Species	Human
GeneID (Human)	7050
GeneID	7050
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-340 of human TGIF1 (NP_733796.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15583
Protein Weight	43kDa
Shipping	Ice bag