


Anti-Tgm1 Rabbit Polyclonal Antibody For Wb, Elisa - P22735

Item Number: CM0018521



Introduction

Rabbit polyclonal antibody against human TGM1 (transglutaminase-1) for WB and ELISA. Detects TGM1 from human, mouse, and rat. Suitable for studies of epidermal cornification, protein cross-linking, and skin barrier function.

[Learn More](#)

Parameter	Value
Product Name	TGM1 Rabbit pAb
Remarks/Alias	LI; KTG; LI1; TKG; ICR2; ARCI1; TGASE; TGM1
Species	Human
Gene ID	7051
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-280 of human TGM1 (NP_000350.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22735
Protein Weight	90kDa
Shipping	Ice bag