

Anti-Tgm3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q08188

Item Number: CM0026877



Introduction

Rabbit polyclonal antibody against human TGM3 (Transglutaminase 3). Suitable for WB, IF/ICC, and ELISA. Detects endogenous TGM3 protein (~77 kDa). Cross-reacts with human and mouse TGM3. Ideal for cornification and skin biology studies.

[Learn More](#)

Parameter	Value
Product Name	TGM3 Rabbit pAb
Remarks/Alias	TGE; UHS2; TGM3
Species	Human
Gene ID (Human)	7053
Gene ID	7053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TGM3 (NP_003236.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08188
Protein Weight	77kDa
Shipping	Ice bag