

Anti-Tgm3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q08188

Item Number: CM0021739



Introduction

Rabbit polyclonal antibody targeting human TGM3 (Transglutaminase 3). Validated for WB, IF/ICC, and ELISA, with cross-reactivity to mouse. Ideal for studies of keratinocyte differentiation, cornified envelope formation, and hair follicle hardening.

[Learn More](#)

Parameter	Value
Product Name	TGM3 Rabbit pAb
Remarks/Alias	TGM3; TGE; UHS2; transglutaminase 3
Species	Human
Gene ID (Human)	7053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 291-480 of human TGM3 (NP_003236.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08188
Protein Weight	76kDa
Shipping	Ice bag