

Anti-Tgm5 Polyclonal Antibody For Wb, Ihc-P, Elisa - O43548

Item Number: CM0026722

Introduction

High-affinity rabbit polyclonal antibody against human TGM5 (Transglutaminase 5), validated for WB, IHC-P, and ELISA. Suitable for cross-linking and cornified envelope research in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	TGM5 Rabbit pAb
Remarks/Alias	TGX; PSS2; TGM6; TGMX; TG(X); TGASE5; TGASEX; TGM5
Species	Human
GeneID (Human)	9333
GeneID	9333
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 461-720 of human TGM5 (NP_963925.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43548
Protein Weight	81kDa
Shipping	Ice bag