


Anti-Tgn46/Tgoln2 Monoclonal Antibody For Wb, Ihc-P, Elisa - O43493

Item Number: CM0003870



Introduction

Rabbit monoclonal antibody against human TGN46 (TGOLN2), suitable for WB, IHC-P, and ELISA. Useful for studying trans-Golgi network membrane trafficking and organelle markers.

[Learn More](#)

Parameter	Value
Product Name	TGN46/TGOLN2 Rabbit mAb
Remarks/Alias	TGN38; TGN46; TGN48; TGN51; TTN2; hTGN46; hTGN48; hTGN51; TGN46/TGOLN2
Species	Human
Gene ID (Human)	10618
Gene ID	10618
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TGN46/TGOLN2 (O43493)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O43493
Protein Weight	46kDa
Shipping	Ice bag