


Anti-Tgn46/Tgoln2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43493

Item Number: CM0019293



Introduction

Rabbit polyclonal antibody against human TGN46/TGOLN2, suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Targets a key trans-Golgi network protein involved in membrane traffic regulation.

[Learn More](#)

Parameter	Value
Product Name	TGN46/TGOLN2 Rabbit pAb
Remarks/Alias	TGN38, TGN46, TGN48, TGN51, TTGN2, hTGN46, hTGN48, hTGN51, TGN46/TGOLN2
Species	Human
Gene ID (Human)	10618
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 323-423 of human TGN46/TGOLN2 (NP_006455.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43493
Protein Weight	46kDa
Shipping	Ice bag