

Anti-Thada Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Yhu6

Item Number: CM0025306



Introduction

High-quality rabbit polyclonal antibody against human THADA (thyroid adenoma-associated protein, GITA/ARMC13). Validated for WB, IF/ICC, ELISA; cross-reacts with human and mouse. Ideal for studying tRNA 2'-O-methylation and thyroid adenoma biology.

[Learn More](#)

Parameter	Value
Product Name	THADA Rabbit pAb
Remarks/Alias	GITA; ARMC13; Trm732; THADA
Species	Human
Gene ID (Human)	63892
Gene ID	63892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1684-1953 of human THADA (NP_071348.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6YHU6
Protein Weight	220kDa
Shipping	Ice bag