

## Anti-Thap2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0W7

Item Number: CM0030449



### Introduction

Rabbit polyclonal antibody against human THAP2, validated for Western blot and ELISA. Cross-reacts with mouse and rat. THAP2 is a zinc-finger DNA-binding protein with roles in transcriptional regulation. Ideal for chromatin research.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | THAP2 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | THAP2                                                                                                                                  |
| Species          | Human                                                                                                                                  |
| GeneID (Human)   | 83591                                                                                                                                  |
| GeneID           | 83591                                                                                                                                  |
| Immunogen        | Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 83-228 of human THAP2 (NP_113623.1) |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, ELISA                                                                                                                              |
| Cross Reactivity | Mouse, Rat                                                                                                                             |
| SWISS            | Q9H0W7                                                                                                                                 |
| Protein Weight   | 26kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |