

Anti-Thap3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wtv1

Item Number: CM0012931

Introduction

Rabbit polyclonal antibody against human THAP3 (Gene ID: 90326), validated for Western blotting and ELISA. Cross-reacts with mouse and rat. Suitable for studying THAP1/THAP3-HCFC1-OGT complex regulation of RRM1 transcription. Weight: 27 kDa.

[Learn More](#)

Parameter	Value
Product Name	THAP3 Rabbit pAb
Remarks/Alias	THAP3
Species	Human
Gene ID (Human)	90326
Gene ID	90326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human THAP3 (NP_001182682.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WTV1
Protein Weight	27kDa
Shipping	Ice bag