

Anti-Thbs1 Polyclonal Antibody For Ihc-P, Elisa - P07996

Item Number: CM0022722



Introduction

Anti-THBS1 rabbit polyclonal antibody validated for IHC-P and ELISA, with cross-reactivity to mouse. Targets 129kDa thrombospondin-1, a multifunctional adhesive glycoprotein involved in angiogenesis and inflammation.

[Learn More](#)

Parameter	Value
Product Name	THBS1 Rabbit pAb
Remarks/Alias	TSP; THBS; TSP1; TSP-1; THBS-1; THBS1
Species	Human
GeneID (Human)	7057
GeneID	7057
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human THBS1 (NP_003237.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P07996
Protein Weight	129kDa
Shipping	Ice bag