

Anti-Thbs1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P07996

Item Number: CM0020535



Introduction

Thrombospondin-1 (THBS1) rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Recombinant immunogen corresponds to amino acids 19-352 of human THBS1. Ideal for studies of angiogenesis, inflammation, and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	THBS1 Rabbit pAb
Remarks/Alias	TSP; THBS; TSP1; TSP-1; THBS-1; THBS1
Species	Human
GeneID (Human)	7057
GeneID	7057
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 19-352 of human THBS1 (NP_003237.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07996
Protein Weight	129kDa
Shipping	Ice bag