

Anti-Thbs1 Rabbit Monoclonal Antibody For Wb - P07996

Item Number: CM0010494



Introduction

High-quality anti-THBS1 rabbit monoclonal antibody validated for WB, IP, IHC-P, and ELISA. Detects human, mouse, and rat THBS1 (UniProt P07996). Ideal for inflammation, angiogenesis, and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	THBS1 Rabbit mAb
Remarks/Alias	TSP; THBS; TSP1; TSP-1; THBS-1
Species	Human
GeneID (Human)	7057
GeneID	7057
Immunogen	A synthetic peptide corresponding to amino acids 19-34 of human THBS1 (NP_003237.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07996
Protein Weight	129 kDa/120 kDa
Shipping	Ice bag