

Anti-Thbs2 Polyclonal Antibody For Wb, Elisa - P35442

Item Number: CM0030377



Introduction

Rabbit polyclonal antibody against human THBS2 (Thrombospondin-2). Suitable for WB and ELISA applications. Cross-reacts with human, mouse, and rat. For cell adhesion and anti-angiogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	THBS2 Rabbit pAb
Remarks / Alias	TSP2, THBS2
Species	Human
Gene ID (Human)	7058
Gene ID	7058
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human THBS2 (NP_003238.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35442
Protein Weight	130 kDa
Shipping	Ice bag