


Anti-Thioredoxin 1 (Trx1/Txn) Monoclonal Antibody For Wb, Ihc-P, Elisa - P10599

Item Number: CM0011640



Introduction

Rabbit monoclonal antibody against human Thioredoxin 1 (Trx1/TXN), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Involved in redox regulation, apoptosis, and transcriptional control.

[Learn More](#)

Parameter	Value
Product Name	Thioredoxin 1 (Trx1/TXN) Rabbit mAb
Remarks/Alias	TRX; TRDX; TRX1; Trx80; Thioredoxin 1 (Trx1/TXN)
Species	Human
GeneID (Human)	7295
GeneID	7295
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 26-105 of human Thioredoxin 1 (Trx1/TXN) (P10599).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10599
Protein Weight	12 kDa
Shipping	Ice bag