


Anti-Thioredoxin 1 (Trx1/Txn) Rabbit Polyclonal Antibody - P10599

Item Number: CM0028863



Introduction

High-quality rabbit polyclonal antibody targeting human Thioredoxin 1 (Trx1/TXN). Validated for WB, IHC-P, IF/ICC, IP, and ELISA applications; cross-reacts with mouse. Suitable for redox signaling and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Thioredoxin 1 (Trx1/TXN) Rabbit pAb
Remarks/Alias	TRX; TRDX; TRX1; Trx80; Thioredoxin 1 (Trx1/TXN)
Species	Human
GeneID (Human)	7295
GeneID	7295
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10599
Protein Weight	12kDa
Shipping	Ice bag