


Anti-Thioredoxin 2 (Txn2) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q99757

Item Number: CM0008902



Introduction

Thioredoxin 2 (TXN2) rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Essential for mitochondrial ROS homeostasis and apoptosis research. Ideal for IVD application development.

[Learn More](#)

Parameter	Value
Product Name	Thioredoxin 2 (Trx2/TXN2) Rabbit mAb
Remarks/Alias	TXN; MTRX; TRX2; MT-TRX; COXPD29; Thioredoxin 2 (Trx2/TXN2)
Species	Human
Gene ID (Human)	25828
Gene ID	25828
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 67-166 of human TXN2 (Q99757)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99757
Protein Weight	18kDa
Shipping	Ice bag