

Anti-Thioredoxin Reductase 2 (Txnrd2) Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nnw7

Item Number: CM0008079

Introduction

Rabbit monoclonal antibody targeting human Thioredoxin reductase 2 (TXNRD2) for WB and ELISA. Detects endogenous TXNRD2 in human, mouse, and rat tissues; important for mitochondrial redox research. SWISS: Q9NNW7.

[Learn More](#)

Parameter	Value
Product Name	Thioredoxin reductase 2 (TXNRD2) Rabbit mAb
Remarks/Alias	TR; TR3; SELZ; GCCD5; TRXR2; TR-BETA; Thioredoxin reductase 2 (TXNRD2)
Species	Human
GenelD (human)	10587
GenelD	10587
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NNW7
Protein Weight	57kDa
Shipping	Ice bag