

Anti-Thiosulfate Sulfurtransferase (Tst) Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q16762

Item Number: CM0017286



Introduction

Rabbit polyclonal antibody targeting human TST (thiosulfate sulfurtransferase/rhodanese), suitable for WB, IHC-P, IP, and ELISA with cross-reactivity to mouse and rat. Ideal for cyanide detoxification and mitochondrial import studies.

[Learn More](#)

Parameter	Value
Product Name	TST Rabbit pAb
Remarks/Alias	RDS, TST
Species	Human
Gene ID (Human)	7263
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human TST (NP_003303.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16762
Protein Weight	33kDa
Shipping	Ice bag