

Anti-Thnsl1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8lyq7

Item Number: CM0021843

Introduction

THNSL1 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA applications. Cross-reacts with human, mouse, and rat. Useful for research on threonine synthase-like 1 in acetylation and proteomics studies.

[Learn More](#)

Parameter	Description
Product Name	THNSL1 Rabbit pAb
Remarks/Alias	TSH1; THNSL1
Species	Human
GeneID	79896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 534-743 of human THNSL1 (NP_079114.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IYQ7
Protein Weight	83kDa
Shipping	Ice bag