


Anti-Thoc3 Polyclonal Antibody For Wb, Elisa - Q96J01

Item Number: CM0013572



Introduction

Rabbit polyclonal antibody targeting THOC3, a component of the TREX complex involved in mRNA export. Validated for Western blot and ELISA, with cross-reactivity to human and mouse. Ideal for studying mRNA processing and nuclear export mechanisms.

[Learn More](#)

Parameter	Value
Product Name	THOC3 Rabbit pAb
Remarks/Alias	THO3; hTREX45; THOC3
Species	Human
GeneID (Human)	84321
GeneID	84321
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96J01
Protein Weight	39kDa
Shipping	Ice bag