


Anti-Thoc4/Allyref Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q86V81

Item Number: CM0030835



Introduction

THOC4/ALYREF Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for mRNA export and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	THOC4/ALYREF Rabbit pAb
Remarks/Alias	ALY; BEF; REF; THOC4; ALY/REF; THOC4/ALYREF
Species	Human
Gene ID	10189
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human THOC4/ALYREF (NP_005773.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86V81
Protein Weight	27kDa
Shipping	Ice bag