


Anti-Thoc6 Polyclonal Antibody For Wb - Q86W42

Item Number: CM0019183



Introduction

Rabbit polyclonal antibody targeting human THOC6 (WDR58/fSAP35) for Western blot and ELISA. Validated in human and mouse samples. Suitable for mRNA processing and nuclear export research.

[Learn More](#)

Parameter	Value
Product Name	THOC6 Rabbit pAb
Remarks/Alias	WDR58, fSAP35, MMRFCGU, THOC6
Species	Human
Gene ID	79228
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human THOC6 (NP_077315.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86W42
Protein Weight	38kDa
Shipping	Ice bag