


Anti-Thoc7 Polyclonal Antibody For Wb, Elisa - Q6I9Y2

Item Number: CM0017650



Introduction

High-quality rabbit polyclonal antibody against THOC7, suitable for Western blot and ELISA. Validated in human, mouse, and rat. Ideal for studying mRNA export and TREX complex function.

[Learn More](#)

Parameter	Value
Product Name	THOC7 Rabbit pAb
Remarks/Alias	fSAP24; hTREX30; NIF3L1BP1; THOC7
Species	Human
Gene ID	80145
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6I9Y2
Protein Weight	24kDa
Shipping	Ice bag